



Teknokroma Capillary Columns

Meta.X5 TRIAZINE (proprietary phase)

Silphenylene phase, selectivity similar to TRB-5, bonded and crosslinked phase.

- New formulation for Meta.X5 stationary phase. Ideal for separation of Triazine herbicides from EPA 609 method.
- Low bleed and excellent inertness for the analysis of traces of herbicides by GC/MS.
- General purpose column for pesticides.

Meta.X5 TRIAZINE

Internal Length	Film	Temp	Part.
Diam.(mm) (m)	Thickness (µm)	limits (°C)	N°. (P/N)
0,30 30	0,25	325 to 350°C	TR-410232

Meta.XLB (proprietary phase)

Silphenylene phase, bonded and crosslinked

- Low polarity phase with Extreme Low Bleed.
- Directly replace for DB-XLB
- General purpose column with extended temperature range (30 to 340/360°C)
- Ideal column for GC-MS analysis
- Unique selectivity for aromatic compounds (PCBs,PAHs,PBDEs)
- Excellent column for pesticides and herbicides

Meta.XLB Equivalent Phase


Restek: Rxi-XLB
Agilent/JW: DB-XLB, VF-XMS
Supelco: MDN 12
Phenomenex: ZB-XLB
Macherey-Nagel: OPTIMA XLB

Triazine Herbicides

Column: **Meta.X5 TRIAZINE**, 30m x 0.25mm x 0.25 µm, P/N TR-410232
 Injection: split 1:25, 250 °C
 Carrier gas: He, ct flow 1.0 ml/min
 Oven: 80 °C (0.5 min) to 160 °C (7 min) @ 30 °C/min to 195°C @ 7 °C/min to 290 °C (3 min) @ 45 °C/min
 Transfer Line: 290°C
 Ionization mode: EI
 Scan range: 50-450 amu
 Sample: Triazine herbicides EPA 619 2 ng/compound on column

Peak Name

- Atraton
- Simazine
- Prometon
- Atrazine
- Propazine
- Terbutylazine
- Secbumeton
- Simetryn
- Ametryn
- Prometryn
- Terbutryn




Meta.XLB

Internal Length	Film	Temp	Part.	
Diam.(mm) (m)	Thickness (µm)	limits (°C)	N°. (P/N)	
0,10 10	0,10	30 to 340/360°C	TR-330141	
0,18 20	0,18	30 to 340/360°C	TR-330984	
	30	0,18	30 to 340/360°C	TR-330934
	60	0,18	30 to 340/360°C	TR-330964
0,25 15	0,10	30 to 340/360°C	TR-330112	
	15	0,25	30 to 340/360°C	TR-330212
	15	1,00	30 to 340/360°C	TR-331012
	30	0,10	30 to 340/360°C	TR-330132
30	0,25	30 to 340/360°C	TR-330232	
30	0,50	30 to 340/360°C	TR-330532	
30	1,00	30 to 340/360°C	TR-331032	
60	0,25	30 to 340/360°C	TR-330262	
0,32 15	0,25	30 to 340/360°C	TR-330213	
	15	1,00	30 to 340/360°C	TR-331013
	30	0,10	30 to 340/360°C	TR-330133
30	0,25	30 to 340/360°C	TR-330233	
30	0,50	30 to 340/360°C	TR-330533	
30	1,00	30 to 340/360°C	TR-331033	
60	0,25	30 to 340/360°C	TR-330263	
0,53 15	1,50	30 to 320/340°C	TR-331515	
	30	1,50	30 to 320/340°C	TR-331535

Analysis of brominated flame retardants (Polybrominated Diphenyl Ethers, PBDEs)

column: **Meta.XLB**, 15 m x 0.18 mm x 0.072 μ m, P/N TR-332414
 injection: 280°C, splitless w/ surge: pulse 40 psi @ 0.30 min, 50ml/min @ 0.4 min
 oven: 100°C (0.5 min) to 320°C (7 min) @ 30°C/min
 carrier gas: Helium, Ct flow 2 ml/min
 liner: Gooseneck splitless (3 mm) w/ Wool
 Detector: MS (SIM)
 ions monitored: 405.9, 460.6, 483.8, 485.8, 563.9, 643.9, 721.6, 799.7
 transfer line temp.: 300°C
 ionization mode: EI
 sample: 1 μ L Mixture BDE-MXF (0.2 ng on-column) & BDE-209 (5 ng on-column) from Wellington Laboratories




Peak Name

- 1 2,4,4'-Tribromodiphenyl ether (BDE-28)
- 2 2,2',4,4'-Tetrabromodiphenyl ether (BDE-47)
- 3 2,3',4,4'-Tetrabromodiphenyl ether (BDE-66)
- 4 2,2',4,4',6-Pentabromodiphenyl ether (BDE-100)
- 5 2,2',4,4',5-Pentabromodiphenyl ether (BDE-99)
- 6 2,2',4,4',6-Pentabromodiphenyl ether (BDE-85)
- 7 2,2',4,4',5,6'-Hexabromodiphenyl ether (BDE-154)
- 8 2,2',4,4',5,5'-Hexabromodiphenyl ether (BDE-153)
- 9 2,2',3,4,4',5',6'-Heptabromodiphenyl ether (BDE-183)
- 10 2,2',3,3',4,5,5',6,6'-Nonabromodiphenyl ether (BDE-209)

TKG 1272

Polychlorinated biphenyl PCB'S

Column: **Meta.XLB**, 30 m x 0.25 mm x 0.25 μ m, P/N TR-330232
 Injection: 1 μ L PCB'S in hexane (200 ppb each compound), 100 ppb
 2,4,5,6-Tetrachloro-m-xylene, 100ppb Decachlorobiphenyl, splitless (0.75 min), 250°C
 Carrier gas: H₂, ct pressure 20 psi (137.8 kPa)
 Linear Velocity: 70 cm/s @ 120°C
 Oven: 120°C (0.5 min) to 260°C (2 min) @ 25°C/min to 325°C (5 min) @ 25 °C/min
 Detector: ECD, 330°C (make up N₂ to ml/min)



Peak Name

- 1 2,Chlorobiphenyl (PCB 1)
- 2 2,4,5,6,-Tetrachloro-m-xylene (SS)
- 3 2,3-Dichlorobiphenyl (PCB5)
- 4 2,2',5-Trichlorobiphenyl (PCB18)
- 5 2,4',5-Trichlorobiphenyl (PCB31)
- 6 2,2',5,5'-Tetrachlorobiphenyl (PCB52)
- 7 2,2',3,5'-Tetrachlorobiphenyl (PCB44)
- 8 2,3',4,4'-Tetrachlorobiphenyl (PCB66)
- 9 2,2',4,5,5'-Pentachlorobiphenyl (PCB101)
- 10 2,2',3,4,5'-Pentachlorobiphenyl (PCB87)
- 11 2,3,3',4,6'-Pentachlorobiphenyl (PCB110)
- 12 2,2',3,5,5',6-Hexachlorobiphenyl (PCB151)
- 13 2,2',4,4',5,5'-Hexachlorobiphenyl (PCB153)
- 14 2,2',3,4,5,5'-Hexachlorobiphenyl (PCB141)
- 15 2,2',3,4,4',5'-Hexachlorobiphenyl (PCB136)
- 16 2,2',3,4,5,5',6-Heptachlorobiphenyl (PCB187)
- 17 2,2',3,4,4',5',6-Heptachlorobiphenyl (PCB183)
- 18 2,2',3,4,4',5,5'-Heptachlorobiphenyl (PCB180)
- 19 2,2',3,3',4,4',5-Heptachlorobiphenyl (PCB170)
- 20 2,2',3,3',4,4',5,5',6-Nonachlorobiphenyl (PCB26)
- 21 Decachlorobiphenyl (SS)

TKG 1265