



Filtering Equipment



Filtering Equipment

- 47 mm filtration apparatus is recommended for filtration of mobile phase and removal of particles from HPLC solvents.
- Manufactured with first quality glass, tube of DURAN glass from Schott.
- The porosity of the filtration plate is of number 3, which means a nominal pore size of 16-40 micrometers.

Reference	Description
TR-F1000	Complete Filtering Equipment 1.000 ml vessel and 250 ml funnel.
TR-F1002	Complete Filtering Equipment 2.000 ml vessel and 250 ml funnel.
TR-F1010	Filtration vessel, frosted glass and 1.000 ml capacity
TR-F1012	Filtration vessel, frosted glass and 2.000 ml capacity
TR-F1022	Fritted glass support, with screw fitting.
TR-F1016	Aluminium plier for the filtering equipment.
TR-F1018	Glass Funnel with 250 ml capacity for the filtering equipment.