




# Teknokroma Capillary Columns

## TRB-35

**35% Diphenyl - 65% Dimethyl polysiloxane, bonded and crosslinked phase.**

- 35% Diphenyl- 65% dimethyl polysiloxane
- Column of intermediate polarity without cyanopropyl groups in its structure
- Excellent confirmation column



Structure of Poly (dimethyldiphenyl) siloxane

### TRB-35 Equivalent Phase

**Agilent:** HP-35, DB-35

**Supelco:** SPB-35

**Restek:** Rtx-35

**SGE:** BPX-35

**Quadrex:** 007-11

## TRB-35

Internal Diam. (mm)	Length (m)	Film Thickness (µm)	Temp limits (°C)	Part. N°. (P/N)
0,25	15	0,15	-20 to 300/320	<b>TR-351312</b>
	15	0,25	-20 to 300/320	<b>TR-350212</b>
	30	0,15	-20 to 300/320	<b>TR-351332</b>
	30	0,25	-20 to 300/320	<b>TR-350232</b>
	60	0,15	-20 to 300/320	<b>TR-351362</b>
	60	0,25	-20 to 300/320	<b>TR-350262</b>
0,32	15	0,15	-20 to 300/320	<b>TR-351313</b>
	15	0,25	-20 to 300/320	<b>TR-350213</b>
	15	0,50	-20 to 290/310	<b>TR-350513</b>
	30	0,15	-20 to 300/320	<b>TR-351333</b>
	30	0,25	-20 to 300/320	<b>TR-350233</b>
	30	0,50	-20 to 290/310	<b>TR-350533</b>
	60	0,15	-20 to 300/320	<b>TR-351363</b>
	60	0,25	-20 to 300/320	<b>TR-350263</b>
	60	0,50	-20 to 290/310	<b>TR-350563</b>
	0,53	15	0,50	-20 to 260/280
15		1,00	-20 to 260/280	<b>TR-351015</b>
30		0,50	-20 to 260/280	<b>TR-350535</b>
30		1,00	-20 to 260/280	<b>TR-351035</b>
60		0,50	-20 to 260/280	<b>TR-350565</b>
60		1,00	-20 to 260/280	<b>TR-351065</b>

